PRESERVATION SOLUTION

Belzer MPS®

(UW Machine Perfusion Solution)

**Container with perfusion and preservation solution
for organs intended for transplantation**

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 **2274**

Container contents:

Solution for perfusion and cold storage of organs for transplantation Belzer MPS® (UW Machine Perfusion Solution)

DESCRIPTION

Composition is: 1000 ml solution contains:

Calcium chloride (dihydrate)	0.068 g	0,5 mmol/l
HEPES (free acid)	2.38 g	10 mmol/l
Potassium phosphate (monobasic)	3.4 g	25 mmol/l
Mannitol (USP)	5.4 g	30 mmol/l
Glucose (anhydrous), beta D (+)	1.80 g	10 mmol/l
Sodium gluconate	17.45 g	80 mmol/l
Magnesium gluconate (anhydrous)	1.13 g	5 mmol/l
Ribose, D (-)	0.75 g	5 mmol/l
Pentafraction (HES)	50.0 g	50 g/l
Glutathione (reduced form)	0.92 g	3 mmol/l
Adenine (free base)	0.68 g	5 mmol/l
Sodium Hydroxide		qs pH: 7.4
Water for injection		qs 1000 ml

Belzer MPS® (UW Machine Perfusion Solution) perfusion solution is a clear to slightly yellow, sterile, non-pyrogenic, non-toxic solution for the in-vitro flushing and temporary continuous perfusion preservation of the explanted organ organs. This solution has an approximate calculated osmolality of 300 mOsm/kg, a potassium concentration of 25 mmol/l sodium 100 mmol/l, and a pH of approximately 7.4 at room temperature.

ACTIONS

Belzer MPS® (UW Machine Perfusion Solution) machine perfusion solution is to be cooled to approximately 5°C (2°C to 8°C) prior to use in a perfusion machine capable of maintaining temperature within the above specified range.

INTENDED USE

Belzer MPS® (UW Machine Perfusion Solution) is intended to be used for the flushing and continuous hypothermic machine perfusion of organs during their storage, transportation, until transplantation into a recipient.

CONTRAINDICATIONS

There are no known contraindications when used as directed.

WARNING

NOT INTENDED FOR DIRECT INJECTION OR INTRAVENOUS INFUSION.

PRECAUTIONS

Belzer MPS® (UW Machine Perfusion Solution) perfusion solution includes constituents (e.g. Pentafraction (HES)), which have caused hypersensitivity reactions in some patients. Physicians should be stay alert to treat possible reactions.

ADVERSE REACTIONS

Needs for dialysis treatment and primary non-functions after kidney transplantation have been reported in clinical tests of the UW machine perfusion solutions. However, those post-transplantation reactions, could not be attributed solely to solutions.

PREPARATION

Remove overwrap prior to use. Check each bag for leaks by squeezing the container firmly. If a leak is found, discard the solution bag. With the overwrap removed, perform a visual inspection of the solution for visible particulate matter. Do not use the solution if obvious particulate matter, precipitates, or contamination are evident in the solution.

After pre-cooling the organ with vascular flushout using Belzer MPS® (UW Machine Perfusion Solution) or other suitable cooled solutions, the graft can be placed in a suitable perfusion apparatus and machine-perfused according to the manufacturer's (or perfusion center's) protocol. The recommended Belzer MPS® (UW Machine Perfusion Solution) solution volume for pulsating organ perfusion is 1000 ml for one kidney. Donor organs must be flushed free of Belzer MPS® (UW Machine Perfusion Solution) prior to transplantation, according to transplantation centers' protocols.

PACKAGING AND STORAGE

EVA bag (PVC free). Belzer MPS® (UW Machine Perfusion Solution) Solution in 10 one-liter (1000 ml) bags per one carton box. Store Belzer MPS® (UW Machine Perfusion Solution) at a temperature between +2°C and +25°C away from direct sunlight. Avoid excessive heat. Do not freeze the solution, and do not use if frozen.

ATTENTION

This medical device is to be used only by the qualified medical staff.

 Date of production	 Caution! See instruction for use.	 2274	The product meets the requirements of the European Medical Devices Directive MDD. The manufacturer is certified by the Notifying Authority No. 2274
 Expiration date	 Do not use if damaged.		The product includes instruction for use.
 Series number	 Sterile (Aseptic filling)	 +25°C	Restriction of acceptable storage temperatures
REF Catalog number	 Single use		Non-pyrogenic
 Keep away from sunlight.	 Medical Device		Carnamedica Sp. z o.o. ul. Olszynki Grochowskiej 21 lok. U6 04-281 Warsaw, Poland www.carnamedica.com mobile: +48 608 336 159

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