PRESERVATION SOLUTION

Belzer UW[®]
Cold Storage Solution
(University of Wisconsin Solution)

**Container with perfusion and preservation solution
for organs intended for transplantation**

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 **2274**

Container contents:

Solution for perfusion and cold storage of organs for transplantation Belzer UW® Cold Storage Solution (University of Wisconsin Solution)

DESCRIPTION

Composition: 1000 ml solution contains:

Pentafraction (HES)	50 g	50 g/l
Lactobionic Acid (as Lactone)	35.83 g	105 mmol/l
Potassium Phosphate monobasic	3.40 g	25 mmol/l
Magnesium Sulfate heptahydrate	1.23 g	5 mmol/l
Raffinose pentahydrate	17.83 g	30 mmol/l
Adenosine	1.34 g	5 mmol/l
Allopurinol	0.136 g	1 mmol/l
Total Glutathione	0.922 g	3 mmol/l
Potassium Hydroxide	5.61 g	100 mmol/l
Sodium Hydroxide/Hydrochloric Acid		qs pH: 7.4
Water for Injection		qs 1000 ml

Belzer UW® Cold Storage Solution (University of Wisconsin Solution) is a clear to light yellow, sterile, non-pyrogenic solution intended for the flushing and the hypothermic storage of organs (kidney, liver, pancreas). Belzer UW® Cold Storage Solution (University of Wisconsin Solution) has an osmolality of approximately 320 mOsm/kg sodium concentration of 29 mmol/l, potassium concentration of 125 mmol/l and a pH of approximately 7.4 at room temperature.

ACTIONS

Belzer UW® Cold Storage Solution (University of Wisconsin Solution) solution, cooled at +2°C to +6°C, is used for flushing of the isolated organ immediately before its removal from the donor and/or immediately after its removal from the donor. The solution is then left in the organ vasculature during hypothermic storage and transportation. Belzer UW® Cold Storage Solution (University of Wisconsin Solution) is not intended to be used with continuous machine perfusion. Administration of Belzer UW® Cold Storage Solution (University of Wisconsin Solution) at the recommended temperature will effectively cool the organ and lower its metabolic requirements.

INTENDED USE

Belzer UW® Cold Storage Solution (University of Wisconsin Solution) is intended for the flushing and the hypothermic storage of kidney, liver, pancreas or other organs at the time of organ removal from the donor, in preparation for storage, transportation and eventual transplantation into the recipient. Storage time of procured organs intended for transplantation is identical to those specified in the instruction manuals for products with identical chemical composition.

CONTRAINDICATIONS

There are no known contraindications when used as directed.

WARNING

NOT INTENDED FOR DIRECT INJECTION OR INTRAVENOUS INFUSION.

PRECAUTIONS

Prior to the reperfusion the donor organ must be flushed with physiological solution to remove Belzer UW® Cold Storage Solution (University of Wisconsin Solution) to prevent occurrence (in the recipient) of potentially serious cardiovascular complications such as hyperkalemic cardiac arrest or bradyarrhythmia. This is necessary because

of the high concentration of potassium in the solution. Belzer UW® Cold Storage Solution (University of Wisconsin Solution) includes components (allopurinol and pentafraction (HES)) which, when administered individually, have caused hypersensitivity reaction in patients.

ADVERSE REACTIONS

Vascular complications (thrombosis, fistulas, primary non-functions, graft rejections, dialysis needs, pancreatitis events after transplantation) have been reported in clinical tests of solutions with the same composition as Belzer UW® Cold Storage Solution (University of Wisconsin Solution). However, those post-transplantation reactions could not be attributed solely to solutions.

PREPARATION

Cool the solution to 2°C to 6°C. Remove overwrap prior to use. Check each bag for leaks by squeezing the container. If a leak is found, discard solution bag. With the overwrap removed, perform a visual inspection of the solution for visible particulate matter. Do not use the solution if obvious particulate matter, precipitates, or contamination are evident in the solution. Use an administration set to flush the organ with Belzer UW® Cold Storage Solution (University of Wisconsin Solution) solution. Prior to connection to the organ, the solution bag should be suspended at a sufficient height to allow for a steady stream of solution and to produce flow rates of at least 30 ml/min during flushing. Flushing should be continued until the organ is uniformly pale and the effluent is clear.

Suggested minimum volumes:

In situ aortic flush:	Ex vivo infusion:	
Adults 2000 ml – 4000 ml Infants 50 ml/kg	<ul style="list-style-type: none"> • Liver Adults 1200 ml Infants 50 ml/kg 	<ul style="list-style-type: none"> • Kidney / Pancreas Adults 300-500 ml Infants 150-250 ml

Additional solution should be dispensed into the container holding the organ. The organ must be totally immersed. Seal the container aseptically. The organ storage container should be maintained within a well-insulated transport container. Ice should be used to surround the organ storage container, but should not be used within the container, where the ice could come into direct contact with the organ. Donor organs must be flushed to remove Belzer UW® Cold Storage Solution (University of Wisconsin Solution) prior to transplantation according to transplantation centers' protocols. Those operations must be done by qualified operators.

PACKAGING AND STORAGE

Bag material: EVA (PVC free). Belzer UW® Cold Storage Solution (University of Wisconsin Solution) solution in 10 one-liter (1.000 ml) bags or in 5 two-liter (2.000 ml) bags per one carton box. Store Belzer UW® Cold Storage Solution (University of Wisconsin Solution) at a temperature between +2°C and +25°C away from direct sunlight. Avoid excessive heat. Do not freeze the solution, and do not use if frozen.

ATTENTION

This medical device is to be used only by the qualified medical staff.

 Date of production	 Caution! See instruction for use.	 The product meets the requirements of the European Medical Devices Directive MDD. The manufacturer is certified by the Notifying Authority No. 2274
 Expiration date	 Do not use if damaged.	 The product includes instruction for use.
 Series number	 Sterile (Aseptic filling)	 Restriction of acceptable storage temperatures
 Catalog number	 Single use	 Non-pyrogenic
 Keep away from sunlight.	 Medical Device	 Carnamedica Sp. z o.o. ul. Olszynki Grochowskiej 21 lok. U6 04-281 Warsaw, Poland www.carnamedica.com mobile: +48 608 336 159

Carnamedica Sp. z o.o.
ul. Olszynki Grochowskiej 21 lok. U6
04-281 Warsaw, Poland
Phone: +48 608336159
Fax: +48 22 3078113
Email: office@carnamedica.com

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Belzer UW[®]
Cold Storage Solution
(University of Wisconsin Solution)

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