# **PRESERVATION SOLUTION**

# **CelProtect®**

**Container with perfusion and preservation solution  
for organs intended for transplantation**

## **INSTRUCTIONS FOR USE**

**EN**

**CarnaMedica Sp. z o.o.**  
ul. Olszynki Grochowskiej 21 lok. U6  
04-281 Warsaw, Poland  
Phone: +48 608336159  
Fax: +48 22 3078113  
Email: [office@carnamedica.com](mailto:office@carnamedica.com)

IfU/CPRT/en/2021/03/04

**CE** 2274

## Container contents:

### Solution for perfusion and cold storage of organs for transplantation CelProtect®

#### QUALITATIVE AND QUANTITATIVE COMPOSITION

Composition: 1.000 ml solution contains:

Glutathione	0.921 g	3 mmol/l
Mannitol	10.930 g	60 mmol/l
Lactobionic acid	28.664 g	80 mmol/l
Glutamic acid	2.942 g	20 mmol/l
Sodium hydroxide	4.000 g	100 mmol/l
Calcium chloride dihydrate	0.037 g	0.25 mmol/l
Potassium chloride	1.118 g	15 mmol/l
Magnesium chloride hexahydrate	2.642 g	13 mmol/l
Histidine	4.650 g	30 mmol/l

#### PHARMACEUTICAL FORM

Solution for organ preservation.

Clear, colourless (or slightly yellow) solution.

pH = 7.3

Osmolality: 320 mOsm/kg

#### CLINICAL PARTICULARS

##### Therapeutic indications

Solution for the preservation of thoracic organs (heart and lung) and abdominal organs (kidney, liver, pancreas) during transplantation: from organ removal from the donor, during storage, transportation, and until transplantation into the recipient.

##### Posology and method of administration

The exact flushing method depends on the centre and on whether or not several organs are to be removed at the same time. Flushing is often carried out in two stages: the 1st flush while the organ is still in situ and the 2nd flush once the organ has been removed. The organ is flushed via a cannula inserted into an artery, at a sufficient pressure to obtain a steady stream of solution in order to produce adequate flushing. In the case of liver transplantation the biliary tree is usually flushed out after explantation, before being placed in the storage and transport container. The volume of flushing solution depends on the organ(s).

Suggested Minimum Volumes:

Heart: Adult: 1-2 litres

Infants: 30ml/kg

– Liver: 6-8 litres

– Lung: 4-6 litres

– Kidney: 4-5 litres

– Pancreas: 4 litres.

– Multi-organ: depending on the organs

**Flushing should be continued until the organ is uniformly pale and the effluent is relatively clear.**

Instructions for use and handling of the product:

- The contents of this bag should not be mixed with other solutions.
- Only intact containers should be used.
- There is no need to filter the solution before use.
- For single use only. Do not reuse.

- A careful visual inspection of the solution for turbidity should be carried out. Only clear and colourless or slightly yellow solutions should be used. If any turbidity, precipitates or contamination are evident, the solution must be discarded.
- Not suitable for continuous machine perfusion. The solution may turn yellow during storage. This does not impair the quality and effectiveness of CelProtect®.

After removal from cold storage (2-8°C), the cooled solution should be used immediately. There have been no adequate clinical studies of the use of CelProtect® for organ transplantation in children.

Storage under cold conditions:

The organ is stored at  $5 \pm 3^\circ\text{C}$  in a sterile container appropriately sized for the organ. The organ must be **completely covered** with the cooled solution. The organ storage container must be aseptically sealed. The container should then be placed in at least a second container. This should be filled with ice, but no ice must enter the organ storage container, where the ice could come into direct contact with the organ. Care must be taken to ensure that the organ is aseptically stored in its container, i.e. the inside of the organ storage container, and its contents must remain sterile.

The additional container(s) must be securely sealed. The organ storage containers must then be maintained within a well-insulated transport container. The organ storage containers should be surrounded with ice. Transport times must be as short as possible.

### Contraindications

Hypersensitivity to the active substances or to any of the excipients.

### Special warnings and precautions for use

- Not suitable for direct injection or for intravenous use in the recipient.
- Exclusively for use for the flushing and cold storage of solid organs.
- Theoretical possibility of residual CelProtect® being released into the general circulation of the recipient upon declamping, if no flushing with physiological solution was applied – it could provoke cardiac arrhythmias or hypotension.

### Interaction with other medicinal products and other forms of interaction

There are no known interactions when used as directed.

### Undesirable effects

- Since CelProtect® is not administered to organ recipients themselves, no solution-specific side effects are to be expected.

### Overdose

Not applicable when used as directed.

## PROPERTIES

### Product

CelProtect® is suitable for heart, lung and abdominal organ transplant preservation.

Administration of the solution at the recommended temperature effectively cools the organ, and its metabolic requirements and associated energy consumption are reduced.

In particular, CelProtect® reduces ischaemic reperfusion damage due to the following properties:

- Prevention of oxidative damage caused by free radicals due mainly to the use of reduced glutathione as the antioxidant.
- Prevention of hypothermia-induced cellular oedema due to the use of membrane-impermeable substances, mannitol and lactobionic acid, which hold water in the extracellular compartment due to their osmotic effect.
- Reduced calcium overload thanks to use solution, e.g. low calcium concentration, moderately high potassium concentrations which is also slightly hyperkalaemic and high sodium and magnesium concentrations comparable to the extracellular environment.
- Regeneration of high-energy substances by supplying them with glutamate, a high-energy substrate, which allows energy to be produced in anaerobic conditions.
- Provision of a buffering capacity due to the use of histidine, which prevents tissue acidosis caused by lactic acid accumulation.

**Pharmacokinetic properties**

Not applicable.

**PHARMACEUTICAL PARTICULARS****List of excipients**

Water for injection

Sodium hydroxide (for pH adjustment).

**Incompatibilities**

In the absence of compatibility studies, this medical device must not be mixed with other medicinal products.

**Special precautions for storage**

Protect against light. Store and transport refrigerated (2°C – 8°C). Do not freeze.

**Nature and contents of container**

CelProtect® is a sterile, ready-to-use solution stored in an EVA type bag. Pack size: 10 bags, 1-liter each or 5 bags, 2-liter each.

**Special precautions for disposal and other handling**

Although the packaging materials are produced under aseptic conditions, the outside of the CelProtect® bag is not sterile. The outside of the bag should be decontaminated if CelProtect® is to be poured out of it. As soon as the outer packaging has been removed, the container should be checked for any leaks by squeezing the bag. If a leak is found, the solution must not be used. The labelled side of the bag should be visible during preparation for use. The spike of a standard cystoscopy infusion set should be inserted into the left port with a twisting motion. The infusion line should be clamped until the start of the infusion. The bag of solution should be placed inside an appropriately sized pressure cuff. The cuff should be inflated to apply sufficient pressure to express the fluid. Alternatively, before infusion, the solution container should be suspended at a sufficient height to guarantee a steady stream of the solution. Prior to transplantation the donor organs must be flushed to remove CelProtect® according to transplantation centers protocols.

**ATTENTION**

This medical device is to be used only by the qualified medical staff.

	Date of production		Caution! See instruction for use.		The product meets the requirements of the European Medical Devices Directive MDD. The manufacturer is certified by the Notifying Authority No. 2274
	Expiration date		Do not use if damaged.		The product includes instruction for use.
	Series number		Sterile (Aseptic filling)		Restriction of acceptable storage temperatures
	Catalog number		Single use		Non-pyrogenic
	Keep away from sunlight.		Medical Device		<b>CarnaMedica Sp. z o.o.</b> ul. Olszynki Grochowskiej 21 lok. U6 04-281 Warsaw, Poland www.carnamedica.com mobile: +48 608 336 159

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